ChromoLympho- β Mix

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English

ChromoLympho-B Proliferation Kit for the *in vitro* stimulation of human B-lymphocytes.

Ref.: EKAMP010MS Q.ty: 5 tubes/5 tests **1 TUBE=1 TEST**
Ref.: EKAMP010M Q.ty: 50 tubes/50 tests **1 TUBE=1 TEST**

Material supplied

Each tube contains ChromoLympho-B Proliferation Factor (synthetic oligodeoxynucleotide, 27bases) and IL-2 Human (Recombinant Human Interleukin-2, source: Yeast) premixed and lyophilized.

Applications:

In vitro stimulation of human B-lymphocytes, in cultured peripheral blood cells or cultured bone marrow cells from patients with chronic lymphocytic leukaemia (CLL) and other B-cell lymphoproliferative syndromes.

Protocol

1. Prepare cells in 9 ml of culture medium (1.11×10^6 cells/ml; blood or bone marrow).
2. Before opening the tube, do a short spin of 20 seconds to collect the pellet down.
3. Dissolve lyophilized powder in 1 ml of culture medium (Chromosome M) and ensure solution is uniform by repeated pipette mixing.
4. Add 1 ml of MIX into 9 ml of culture cell preparation (10^6 cells/ml final cell concentration; blood or bone marrow).

Protocol: Synchronized culture

1. After 48-72 H of culture, add 0.2 ml of Synchroset solution A (Ref. EKAMTS-008) into the 10 ml cell culture.
2. After overnight incubation, add 0.2 ml of Synchroset solution B (Ref. EKAMTS-008), without changing the culture medium, and incubate further 5 hours.
3. Add few drops of Colcemid (10 µg/ml- ECM0040B), and incubate at 37 -/+ 2°C during 60 minutes (if Synchroset is used for the first time, different incubation times with Colcemid should be tried, i.e. 30 min, 45 min, 1h).
4. Follow the standard protocol for metaphase preparation.

Storage/Stability:

ChromoLympho-B MIX is shipped at room temperature, and must be stored at -20°C upon arrival. Lyophilized powder is stable at -20°C until expiration date provided on the tube. After reconstitution, store liquid working aliquots at -20°C until use. ChromoLympho-B MIX solution is stable during 6 months at -20°C upon reconstitution in sterile conditions.

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EuroClone®
serving science through innovation

BIOAIR®

EuroClone®
serving science through innovation

EuroClone S.p.A.

Via Figino, 20/22 - 20016 Pero (MI) Italy
010 +39 02 381951 - 02 +39 02 38101465
info@euroclone.it - www.euroclone.it