ChromoLympho-B

Mix

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Protocol of use

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ChromoLympho-B Proliferation Kit for the *in vitro* stimulation of human B-lymphocytes.

Ref.: EKAMP010M Q.ty: 50 tubes/50 tests **1 TUBE=1 TEST**

Material supplied

50 Tubes: each tube contains ChromoLympho-B Proliferation Factor (synthetic oligodeoxynucleotide, 27bases) and IL-2 Human (Recombinant Human Interleukin-2, source: Yeast) premixed and lyophilised.

Applications:

In vitro stimulation of human B-lymphocytes, in cultured peripheral blood cells or cultured bone marrow cells from patients with chronic lymphocytic leukaemia (CLL) and other B-cell lymphoproliferative syndromes.

English

Protocol

1. Prepare cells in 9 ml of culture medium (1.11x10⁶ cells/ml; blood or bone marrow).
2. Before opening the tube, do a short spin of 20 seconds to collect the pellet down.
3. Dissolve lyophilized powder in 1 ml of culture medium (Chromosome M) and ensure solution is uniform by repeated pipette mixing.
4. Add 1 ml of MIX into 9 ml of culture cell preparation (10⁶ cells/ml final cell concentration; blood or bone marrow).

Protocol: Synchronized culture

1. After 48-72 H of culture, add 0.2 ml of Synchronset solution A (Ref. EKAMTS-008) into the 10 ml cell culture.
2. After overnight incubation, add 0.2 ml of Synchronset solution B (Ref. EKAMTS-008), without changing the culture medium, and incubate further 5 hours.
3. Add few drops of Colcemid (10 µg/ml- ECM0040B), and incubate at 37 -/+ 2°C during 60 minutes (if Synchronset is used for the first time, different incubation times with Colcemid should be tried, i.e. 30 min, 45 min, 1h).
4. Follow the standard protocol for metaphase preparation.

Storage/Stability:

ChromoLympho-B MIX is shipped at room temperature, and must be stored at -20°C upon arrival. Lyophilized powder is stable at -20°C until expiration date provided on the tube. After reconstitution, store liquid working aliquots at -20°C until use. ChromoLympho-B MIX solution is stable during 6 months at -20°C upon reconstitution in sterile conditions. Please avoid freeze-thaw cycles.

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EuroClone
serving science through innovation

EuroClone
serving science through innovation

EuroClone S.p.A.

Via Figino 20/22 - 20016 Pero (MI) Italy

☎ + 39 02 38195.1 📠 + 39 02 38101465

✉ info@euroclone.it 🌐 www.euroclone.it