

ENGLISH

For Professional Use Only

NTRK2 Break Apart FISH Probe Kit

Introduction


The NTRK2 Break Apart FISH Probe Kit is designed to detect rearrangements in the human NTRK2 gene located on chromosome band 9g21.33. In addition to revealing breaks, which can lead to translocation of parts of the gene, inversion, or its fusion to other genes, the probe set can also be used to identify other NTRK2 aberrations such as deletions or amplifications. Rearrangements and abnormal expression of the NTRK2 gene - also known as OBHD, TRKB, DEE58, trk-B, EIEE58 or GP145-TrkB - has been observed in neuroblastoma, pancreatic ductal adenocarcinoma, Wilms' tumors and colorectal cancer, and a number of developmental and metabolic disorders.

Intended Use

To detect rearrangements in the human NTRK2 locus situated on chromosome band 9g21.33.

Cont.	Color
LSP NTRK2 5' FISH Probe	CytoOrange
LSP NTRK2 3' FISH Probe	CytoGreen

Probe Design



LSP NTRK2 5' FISH Probe covers the 5' (start) portion of the NTRK2 locus and some adjacent genomic sequences. LSP NTRK2 3' FISH Probe covers sequences at the 3' (end) of the gene. The two probes are flanking sequences across the NTRK2 locus in which variable breakpoints have been observed.

Not to Scale

Cat. No.	Volume
CT-PAC396-10-OG	10 Tests (100 μL)

Signal Pattern Interpretation

Normal Patterns **Abnormal Patterns** 2F* Other Patterns

*Overlapping orange and green signals can appear as yellow.

IVD



CytoTest Inc. 1395 Piccard Drive, Suite 308 Rockville, MD 20850, USA

¹⁾ Nakagawara, A, et al. *Genomics* 25(2):538-46 (1995).
2) Lange, AM & Lo, HW. *Cancers (Basel)* 10(4) pii: E105. doi: 10.3390/cancers10040105 (2018).
3) Prabhakaran, N, et al. *Neuropathology* Mar 4. doi: 10.1111/neup.12458. (2018).
4) Suurmeijer, AJH, et al. *Genes Chromosomes Cancer* 57(12):611-621 (2018).
5) Remoué, A, et al. *Pathol. Int.* 2:94-96 (2019).

^{*} CE IVD only available in certain countries. All other countries are either ASR or RUO. Please contact your local dealer or our headquarters for more information.