

User Guide



Version 1.0 (April 2007)

Product: FISHBright™ Nucleic Acid Labeling Kits

Catalog number (s): FLK-001 – FLK-007

Lot number: See label on package

Unit Size: 10 labeling reactions

For the full manual please visit www.kreatech.com

I. Components & Storage - IMPORTANT

Materials Provided	Quantity	Storage Conditions
ULS hapten/dye	20 µl	2-6°C
10 x Labeling solution	100 µl	2-6°C
KREApure™	10 pcs	2-6°C
Cell Hybridization Buffer	1 ml	2-6°C
Paraffin Hybridization Buffer	1 ml	2-6°C
KREAbboost	100 µl	-20°C*

* KREAbboost is shipped at 4°C, store at -20°C upon delivery

Reagents and Buffers not Included in Kit

DNase I, Nicking buffer, Nuclease free H₂O, NH₄AC, EtOH

II. Suggested Fragmentation Protocols

- If DNA fragmentation is performed using DNase I, it must be done **before** ULS labeling. If fragmenting DNA by sonication then this may be carried out before or after ULS labeling
- For in situ hybridization the optimal probe fragment size is between 100-300 bp. Probes over 1 kb will generate a spotted background in *in situ* hybridization. Because most PCR generated probes are smaller than 1 kb, these probes can be labeled without fragmentation. Other probes such as plasmids, cosmids, YACs, PACs and BACs must be fragmented **prior** to labeling with ULS.

a. Suggested Fragmentation by Sonication

- Prepare at least 100 µL of DNA solution in nuclease free water with a minimal DNA concentration of 50 ng/µL
- Sonicate for 5 cycles of 1 minute each at amplitude 5µm, keep sample on ice. Cool the DNA solution on ice for 1 minute after each cycle. Check our website (www.kreatech.com) for more protocols.
- Check the fragment length by electrophoresis on a 1% agarose gel
- Label the DNA according to the labeling protocol

b. Suggested Fragmentation by DNase I Treatment

- Prepare DNase I stock solution: dissolve 1mg of DNase I (Roche # 104 159, approx. 2000 units/mg) in 1 mL of 5 mM NaAc pH 5.2, 1 mM CaCl₂, 50 % glycerol. Make sure buffer is ice cold before adding the lyophilized DNase I. Mix by inversion. Do not vortex. Store stock DNase I solution at -20°C; avoid freeze-thaw cycles

- Prepare 10 x nicking buffer: 50 mM Tris-HCl pH 7.5, 10 mM CaCl₂, 10 mM MgCl₂
- Dilute the DNase I stock 1:5000 in 1x nicking buffer. The buffer solution should be cold before the addition of the enzyme. These dilutions should be prepared immediately before use
- Setup nicking reaction on ice as follows:
 - 1 µg template DNA
 - 2.5 µL 10 x nicking buffer
 - 3-5 µL diluted DNase I
 - H₂O to 25 µL
- Incubate at 37°C for 10 minutes. Stop reaction by placing on ice for at least 1 minute and precipitate with ¼ volume of 10M NH₄Ac and 2.5 volumes of 100% EtOH. Resuspend the Nuclease-free water
- Check the fragment length by electrophoresis on a 1% agarose gel

III. ULS Labeling

For labeling 1µg nucleic acid sample

- Briefly spin all required reagents to collect contents of tubes
- Take 1 µg of nucleic acids sample **ensure final concentration in labeling reaction is above 50 ng/µL**. (If the final concentration of the nucleic acid sample in the labeling mixture is below 50 ng/µL then the optimal degree of modification is not achieved)
- Add 2µL of ULS label per 1µg nucleic acid sample
- Add 1/10 volume of 10x Labeling solution
- Adjust with Nuclease-free Water to final volume and mix by pipetting (see example set-up below)
- Label sample by incubation for 30 minutes at 85°C
- Place samples on ice, spin down to collect contents of tube before proceeding with purification using the KREApure columns

Example of ULS labeling of 1 µg nucleic acid sample in 20µL

Nucleic acid sample in Nuclease-free Water	16 µL
ULS dye	2 µL
10 x labeling solution	2 µL
Total volume	20 µL

Dye removal using KREApure columns

(20800 x g is equivalent to 14,000 rpm on eppendorf 5417C)

- Resuspend column material by vortexing
- Loosen cap ¼ turn and snap off the bottom closure
- Place the column in a 2 mL collection tube
- Pre-spin the column for 1 minute at 20800 x g
- Discard flow through and re-use collection tube
- Wash the column with 300 µL nuclease-free water
- Discard collection tube and flow-through

