# *Amniomed* *plus*

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**Protocol of use**

**Protocolo de utilización**

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# English

## Complete medium for culture of amniocytes and chorionic villi

### Material supplied

EK AMG 200, 1 bottle of 100 ml  
EK AMG 600, 1 bottle of 500 ml

### Storage and stability

Amniomed Plus must be stored at  $-20 \pm 2^\circ\text{C}$ . Thaw the Amniomed Plus in a  $37^\circ\text{C}$  water bath and mix well by swirling prior to use. Addition of L-Glutamine, antibiotics and serum are not necessary since these components are already present. Once thawed, it must be stored at a temperature between  $+2^\circ\text{C}$  and  $+8^\circ\text{C}$  for maximum 8 days. The shelf life of the product is 2 years.

### Protocol of use

NB: In situ culture protocols refer to the use of Amnioslide. With regards to the use of Amnioidish refer to the values in brackets.

### AMNIOCYTES - Sample treatment

1. Transfer the amniotic fluid into sterile tubes and centrifuge 10 minutes at 1000 rpm.
2. Remove the supernatant and resuspend the pellet in 2-3 ml of Amniomed Plus.

#### In situ culture

1. Place 1 ml [0,5 ml] of cell suspension gently onto each Amnioslide, on the microscope slide. Incubate a  $37^\circ\text{C}$ , 5%  $\text{CO}_2$  for 6- 12 hours, with minimal mechanical disturbance.
2. Add 4 ml [2 ml] of fresh Amniomed Plus.
3. On Day 6, replace the culture medium with 5 ml [2,5 ml] of fresh Amniomed Plus.
4. On Day 7/Day 8 check for progress of growth using an inverted microscope. At the optimal cell growth, add 30  $\mu\text{l}$  [15  $\mu\text{l}$ ] of Colcemid (10  $\mu\text{g/ml}$ ).
5. Incubate at  $37^\circ\text{C}$ , 5%  $\text{CO}_2$  for 60 - 90 minutes.

#### Flask culture

1. Prepare 2 to 4 sterile flasks of  $25\text{ cm}^2$ . Add 4 ml of Amniomed Plus to each flask.
2. Add 0.5 ml of suspension and incubate at  $37^\circ\text{C}$ , 5%  $\text{CO}_2$ .
3. On the 6th day, completely change the culture medium with 5 ml of Amniomed Plus.
4. On the 7th/8th day or when suitable growth has been reached, add 30  $\mu\text{l}$  Colcemid (10  $\mu\text{g/ml}$ ).
5. Incubate at  $37^\circ\text{C}$ , 5%  $\text{CO}_2$  for 60 – 90 minutes.

### CHORIONIC VILLI (CVS) - Sample treatment

1. Transfer the specimen into a 60 mm Petri dish containing 5 ml of RPMI 1640.
2. Wash villi with fresh medium to remove blood cells.
3. Using the inverted microscope, carefully dissect any remaining clots or decidual fragments.
4. Transfer the villi in a 15 ml sterile centrifuge tube containing 1 ml of Pronase E ( $4 \times 10^6$  PU/g) and incubate at room temperature for 4-6 minutes.
5. Add 3-5 ml of cold Hank's Balanced Salts Solution (kept at  $+4^\circ\text{C}$ ).
6. Centrifuge for 5 minutes at 1500 rpm and discard the supernatant.
7. Add 2 ml of sterile Collagenase type II (1 mg/ml) and incubate at  $37^\circ\text{C}$  for 10 minutes.
8. Add 3-5 ml of cold HBSS (kept at  $+4^\circ\text{C}$ ).
9. Centrifuge for 5 minutes at 1500 rpm and discard the supernatant.
10. Add 2 ml of Amniomed Plus and resuspend the pellet.

#### In situ culture

1. Prepare 2 to 6 Amnioslide, depending on the number of cells. Add 4,5 ml [2 ml] of Amniomed Plus medium to each Amnioslide.
2. Drop about 500  $\mu\text{l}$  of the cell suspension on each Amnioslide.
3. Incubate at  $37^\circ\text{C}$  in 5%  $\text{CO}_2$ .
4. On Day 5 replace the culture medium with 5 ml [2,5 ml] of fresh Amniomed Plus.
5. On Day 6 / Day 7 check for progress of growth using an inverted microscope. At the optimal cell growth, add 30  $\mu\text{l}$  [15  $\mu\text{l}$ ] of Colcemid (10  $\mu\text{g/ml}$ ).
6. Incubate at  $37^\circ\text{C}$ , 5%  $\text{CO}_2$  for 4 - 6 hours.

#### Flask culture

See protocol for the amniocytes culture in flask. Change the timing of points 3 and 4 according to the in situ CVS culture protocol, points 3 and 4.

### Chromosomes fixation and slide preparation from in situ culture

1. Suck the medium up completely.
2. Drop 6 ml [3 ml] of hypotonic solution (distilled  $\text{H}_2\text{O}$ , 0,6% sodium citrate, 0,1% KCl) and incubate at room temperature for 10 minutes.
3. Suck the hypotonic solution up and drop 6 ml [3 ml] of Ibraimov solution (distilled  $\text{H}_2\text{O}$ , 5% acetic acid), incubate at room temperature for 5 minutes.
4. Suck the Ibraimov modified solution up and drop 6 ml [3 ml] of fresh fixative mixture (methanol:acetic acid, 4:1) at room temperature. The time for this wash has no influence on final results.
5. Repeat 4) twice.
6. Remove the microscope slide and dry under constant condition of right temperature and relative humidity into Optichrome (ref. EKAMH950).

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**EuroClone S.p.A.**
EuroClone S.p.A. has a Quality System certified in compliance