



EuroClone®

serving science through innovation

Composition of Trypsin 2.5 % in HBSS w/o Phenol Red w/o Calcium w/o Magnesium

Code: ECB3051

Inorganic salt :

Potassium chloride	400 mg/l
Potassium phosphate monobasic anhydre	60 mg/l
Sodium bicarbonate	350 mg/l
Sodium chloride	8000 mg/l
Sodium phosphate dibasic anhydrous	47.88 mg/l

Other components :

D-Glucose	1000 mg/l
Trypsin 1:250	25000 mg/l

EuroClone S.p.A.

Via Figino 20/22, 20016 Pero (MI) - Italy
Phone: +39.02.38195.1 - Fax +39.02.38101465
e-mail: info@euroclone.it



IOAIR®



www.euroclone.it